


MALE URINAL

H140-01

*Supports your initiative to reduce
Urinary Tract Infections (UTI)*

FEATURES/BENEFITS:

- Improved smooth rolled-edge opening
- Foot design for stability empty or full
- Translucent urinals for content visualization and measuring ease
- Specially designed handle conveniently hangs on most beds or can be used with the Uri Clip™ H143-01
- 1000cc/32 oz capacity with clear graduations
- Reduction in packaging by 12% = Improved customer stocking and warehousing efficiencies
- Improved manufacturing process, reduction in plastic by 18% = less plastic waste in landfills
- Bulk Packaging 50/cs
- Recyclable 
- **Made in the U.S.A.**



 **MEDEGEN**
Medical Products
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298 • Customer Service Fax 631-714-6311

E-mail cservice@medegenmed.com



To view our products online, please visit us at www.medegenmed.com or directly with your device by scanning the QR code to the left.



FEATURES/ BENEFITS:

- Single Patient Use and/or Reusable
- Sterilizable options
- Translucent urinals for content visualization and measurement ease
- Wide base design for stability empty or full
- Odor shield (lid) options

Packaging Options:

- Bulk
- Individually bagged

Medegen Medical Products

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298

E-mail cservice@medegenmed.com

Web www.medegenmed.com

Urinals - Male (Single Patient Use)

Code#	Color	Capacity	Qty/cs	Packaging
H140-01	Translucent	1qt/1000	50	Bulk
H140CR-01	Translucent w/Comfort Ring	1qt/1000	50	Bulk
H140D-01	Translucent	1qt/1000	50	Bulk
H141-01	Translucent	1qt/1000	50	Individually Bagged

Urinals - Male (Sterilizable*)

Code#	Color	Capacity	Qty/cs	Packaging
00090	Blue	1.25qt/1180cc	12	Bulk w/Liner
00095	Blue	1qt/1000cc	12	Individually Bagged
00195	Translucent w/Blue Cover	1qt/1000cc	12	Individually Bagged
89150	Stainless Steel		6	

*Sterilizable products can withstand temperature up to 275°F (135°C), steam sterilization

Urinals - Female (Single Patient Use)

Code#	Color	Capacity	Qty/cs	Packaging
H145-01	Translucent	1qt/1000cc	36	Bulk

Uri-Clip™ Urinal Holder

Code#	Color	Capacity	Qty/cs	Packaging
H143-01	White	N/A	48	Bulk