

FEMALE URINAL


H145-01

*Supports your initiative to reduce
Urinary Tract Infections (UTI)*

A Better Alternative to a Fracture Bedpan!

- More comfortable for the patient
- Patient can use at their own convenience
- Nurse time saver
- Less prone to spillage
- Comfortable, elongated body design

FEATURES/BENEFITS:

- Contoured smooth opening
- Contents easily measured
- Specially designed handle conveniently hangs on most beds or with H143-01 Uri-Clip™
- Pronounced neck angle helps confine contents when urinal is positioned horizontally
- Graduated: 25 cc and 1 oz. increments up to 1000 cc/32 oz
- Packaging: 36/case, bulk
- Recyclable 



 **MEDEGEN**
Medical Products
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036
Toll Free 800-511-6298 • Customer Service Fax 631-714-6311
E-mail cservice@medegenmed.com



To view our products online, please visit us at www.medegenmed.com or directly with your device by scanning the QR code to the left.



FEATURES/ BENEFITS:

- Single Patient Use and/or Reusable
- Sterilizable options
- Translucent urinals for content visualization and measurement ease
- Wide base design for stability empty or full
- Odor shield (lid) options

Packaging Options:

- Bulk
- Individually bagged

Medegen Medical Products

209 Medegen Drive, Gallaway, TN 38036

Toll Free 800-511-6298

E-mail cservice@medegenmed.com

Web www.medegenmed.com

Urinals – Male (Single Patient Use)

Code#	Color	Capacity	Qty/cs	Packaging
H140-01	Translucent	1qt/1000	50	Bulk
H140CR-01	Translucent w/Comfort Ring	1qt/1000	50	Bulk
H140D-01	Translucent	1qt/1000	50	Bulk
H141-01	Translucent	1qt/1000	50	Individually Bagged
H149-01	Translucent	1qt/1000	50	Bulk

Urinals – Male (Sterilizable*)

Code#	Color	Capacity	Qty/cs	Packaging
00090	Blue	1.25qt/1180cc	12	Bulk w/Liner
00095	Blue	1qt/1000cc	12	Individually Bagged
00195	Translucent w/Blue Cover	1qt/1000cc	12	Individually Bagged
89150	Stainless Steel		6	

*Sterilizable products can withstand temperature up to 275°F (135°C), steam sterilization

Urinals – Female (Single Patient Use)

Code#	Color	Capacity	Qty/cs	Packaging
H145-01	Translucent	1qt/1000cc	36	Bulk

Uri-Clip™ Urinal Holder

Code#	Color	Capacity	Qty/cs	Packaging
H143-01	White	N/A	48	Bulk