



# AUTOCLAVE BAGS

*Autoclavable Biohazardous Waste Bags withstand the rigors of high temperature / high pressure steam while meeting requirements for storage, transport, and processing of biohazardous waste within the facility.*

- **Dual-Tested Autoclave Bags** are ASTM D1922 tested (Elmendorf Tear) to exceed 480g of tear strength (MD, TD) and ASTM D1709 tested to exceed 165g of dart impact. Complies with US DOT 49CFR 173-197 and regulations in all 50 states. **Item listing at right.**
- **Standard Autoclave Bags** are ASTM D1709 tested to exceed 165g of dart impact. Complies with regulations within states that do NOT have dual-test bag requirements. **Item listing on back.**

## FEATURES/BENEFITS:

- **Operating Temperature** – Range of 250-274°F (121-134°C) with a max of up to 285° F (140°C).
- **Temperature Indicator** – Printed temperature block turns brown when autoclave temperature is reached (250°F, 10 min+/- 5 min).
- **Peace of Mind Performance** – Dual-tested bags are compliant with regulations in all 50 states.
- **Biohazard Symbol** – With wording in easy-to-see contrasting background color.
- **Languages** – Printed in English and Spanish OR French
- **Colors** – Available in red, clear, orange and buff for easy waste stream or department identification.
- **Twist Ties** – Included for fast closure and opening as needed.
- **Custom Printing Available (facility name, etc.)** – contact your local sales representative.



To view our products online, please visit us at [www.medegenmed.com](http://www.medegenmed.com) or directly with your device by scanning the QR code to the left.



Item No.	Size (in)	Gauge	Film Color	Print	Bags/Cs
AC1419C	14x19	2 mil	Clear	Black	500
AC1419R	14x19	2 mil	Red	Black	500
AC1923C	19x23	2 mil	Clear	Black	200
AC1923R	19x23	2 mil	Red	Black	200
AC2530C	25x30	2 mil	Clear	Black	200
AC2530R	25x30	2 mil	Red	Black	200
AC2535C	25x35	2 mil	Clear	Black	200
AC2535R	25x35	2 mil	Red	Black	200
AC3138C	31x38	2 mil	Clear	Black	100
AC3138R	31x38	2 mil	Red	Black	100
AC3847C	38x47	2 mil	Clear	Black	100
AC3847R	38x47	2 mil	Red	Black	100



209 Medegen Drive, Gallaway, TN 38036  
 Toll Free 800-511-6298 • Customer Service Fax 631-714-6311  
 E-mail [cservice@medegenmed.com](mailto:cservice@medegenmed.com)

## Dual-Tested Autoclavable Biohazardous Waste Bags

with Indicator, Meet A.S.T.M. Dart Impact/Elmendorf Tear & DOT CFR 49 173-197 Requirements - Flat Seal, Coreless Roll, Twist Ties, Max. Temp 285°F Print/Label: Biohazard Labeling, Biohazard Symbol, English/Spanish  
SEE PRODUCT LISTING ON FRONT OF SHEET.



## Autoclavable Biohazardous Waste Bags

Polypropylene Film with Indicator, Meet A.S.T.M. Dart Impact Requirements - Flat Seal, Coreless Roll w/Twist Ties, Max. Temp 285°F Print/Label: Biohazard Labeling, Biohazard Symbol, English/Spanish

Code#	Bag/Print Color	Size (inches)	Gauge/Mil	Capacity/gal	Qty/cs	Qty/roll
880	Orange/Black	8 x 12	1.80	1-2	400	20
5001.4	Red/Black	8 x 12	2.0	1-2	400	20
852	Clear/Red	12 x 24	1.80	2-3	200	20
5010	Red/Black	12 x 24	2.0	2-3	500	20
881	Orange/Black	14 x 18	1.80	3-4	400	20
5012	Red/Black	14 x 19	2.0	3-4	500	20
5119	Buff/Red	19 x 23	1.75	7-10	200	20
8-702	Orange/Black	19 x 23	2.0	7-10	200	20
882	Orange/Black	19 x 24	1.80	7-10	200	20
842	Red/Black	19 x 24	1.80	7-10	200	20
5013	Red/Black	19 x 24	2.0	7-10	200	20
8-150	Clear/Red	24 x 30	1.80	12-14	100	10
8-504	White/Black	24 x 33	2.0	12-14	200	20
8-250	Clear/Red	24 x 36	1.80	13-16	100	10
854	Clear/Red	25 x 30	1.80	12-14	200	20
844	Red/Black	25 x 30	1.80	12-14	200	20
846	Red/Black	25 x 35	1.80	13-16	200	20
884	Orange/Black	25 x 35	1.80	13-16	200	20
8-200	Clear/Red	26 x 30	1.80	12-14	100	10
855	Clear/Red	31 x 38	1.80	20-30	100	10
856	Red/Black	31 x 38	1.80	20-30	100	10
858	Buff/Red	31 x 38	1.80	20-30	100	10
5034	Clear/Red	33 x 40	2.0	33	200	20
8-300	Clear/Red	36 x 35	1.80	33	100	10
8-640	Clear/Red	36 x 45	1.80	44	100	10
5055.HTS	Red/Black	38 x 46	1.90	44	100	10
5138	Buff/Red	38 x 46	2.0	44	100	10
5053	Clear/Red	38 x 47	1.75	44	100	10
872	Clear/Red	38 x 47	1.80	44	100	10
848	Red/Black	38 x 47	1.80	44	100	10
886	Orange/Black	38 x 47	1.80	44	100	10
864	Buff/Red	38 x 47	1.80	44	100	10
5028	Red/Black	40 x 70	2.0	55-60	100	10
2120	Clear/Red	45 x 58	1.80	55	100	10
5192	Clear/Red	45 x 58	2.0	55	100	10

## Max. Temp 285°F Print/label: Biohazard Labeling, Biohazard Symbol, English/French

Code#	Bag/Print Color	Size (inches)	Gauge/Mil	Capacity/gal	Qty/cs	Qty/roll
5020.F	Clear/Red	12 x 24	1.75	2-30	500	20
842Y	Yellow/Black	19 x 24	2.20	7-10	200	20
5062.F	Yellow/Black	20 x 22	1.75	7	250	25
5142.F	Clear/Red	24 x 36	1.75	12-16	200	20
5025.F	Clear/Red	25 x 30	1.75	12-14	200	20
844Y	Yellow/Black	25 x 30	2.20	12-14	200	20
846Y	Yellow/Black	25 x 35	2.20	13-16	200	20
5064.F	Yellow/Black	30 x 38	2.0	20-30	125	25
856Y	Yellow/Black	31 x 38	2.20	20-30	100	10
848Y	Yellow/Black	38 x 47	2.20	44	100	10

## No Print/Label

Code#	Bag/Print Color	Size (inches)	Gauge/Mil	Capacity/gal	Qty/cs	Qty/roll
5023	Clear	38 x 47	1.80	44	100	10